

| | | | |
|----------|--|----------|--|
| ANALYST: | | VPDES NO | |
|----------|--|----------|--|

Parameter: Total Kjeldahl Nitrogen
Method: Direct Potentiometric
 04/01

METHOD OF ANALYSIS:

| | |
|--|---|
| | EPA Methods for Chemical Analysis 351.4 |
|--|---|

| | Y | N |
|---|---|---|
| 1) Is analysis performed with a pH meter with an expanded millivolt scale capable of 0.1 mV resolution between - 700 mV and + 700 mV or an ion specific meter? [5.1] | | |
| 2) Is a Teflon coated stirrer bar and magnetic stirrer used during the analysis? [5.3] | | |
| 3) Is a thermally insulated magnetic stirrer used during the analysis? [5.3] | | |
| 4) Is the electrode used an Orion Model 95-12 or 95-10, EIL Model 8002-2, Beckman Model 39565 or equivalent? [5.2] | | |
| 5) For short term storage (week or less) is probe stored as specified in the manufacturers instructions? [Mfr.] | | |
| 6) For long term storage (longer than one week) is the electrode drained, rinsed with distilled water and stored dry? [Mfr.] | | |
| 7) Is ammonia free water prepared immediately before use and used for all aspects of the procedure? [6.1] | | |
| 8) Are standards prepared using Class A volumetric glassware? [Permit] | | |
| 9) Are standards treated the same as samples? [7.4.1] | | |
| 10) For macro Kjeldahl analysis is 15 mL of 10N NaOH added to 100 mL of digested sample (micro Kjeldahl 6 mL 10N NaOH to a 50 mL aliquot)? [macro-7.4.2; micro-7.2.1 and 4.4.3] | | |
| 11) Is 4 mL of NaOH-Nal-EDTA reagent added to the sample after the electrode has been immersed? [7.4.2 / 7.4.3] | | |
| 12) Is the calibration curve prepared using 4 cycle semilogarithmic paper, with the ammonia nitrogen in mg/L on the logarithmic axis and the electrode potential in mV on the linear scale; or by using a computer program which has been verified by either a hand-held calculator or semilog paper? [8] | | |
| 13) Are direct readout meters calibrated according to manufacturer's instructions? [Mfr.] | | |
| 14) Are standards read from lowest to highest concentration? [8] | | |
| 15) Is the instrument slope documented to be within manufacturer's specifications each sample run? (Corning -55 \pm 5 mV, Orion -57 \pm 3 mV, Accumet -59 \pm 4 mV, Hach -58 \pm 4 mV) [Mfr.] | | |
| 16) Is the electrode rinsed with distilled water and blotted dry between measurements? [Mfr.] | | |
| 17) Are samples and standards stirred so that bubbles are not sucked in the solution? [Mfr.] | | |
| 18) Is the electrode held at a 20 - 30 degree angle in the sample during analysis? [Mfr.] | | |
| 19) Is the electrode tip free of bubbles during operation? [Mfr.] | | |
| 20) Is a new curve drawn when calibration standards are not within \pm 5.0% of the curve? [Permit] | | |

| | Y | N |
|--|---|---|
| 21) Are results recorded as soon as the meter stabilizes? [7.4] | | |
| 22) Are results recorded in terms of ammonia nitrogen? [8] | | |
| 23) Are standards and samples at the same temperature when analyzed? [7.4.1] | | |
| 24) Are all calibrations, calculations, temperatures recorded? [Permit] | | |

PROBLEMS: