

<b>ANALYST</b>		<b>LAB CODE</b>	
----------------	--	-----------------	--

Parameter: Base/Neutrals and Acids (GC/MS)

Checklists

04/02

INSTRUMENT \_\_\_\_\_

DATE \_\_\_\_\_

**METHOD OF ANALYSIS**

	EPA 625 Series, 40 CFR, Part 136
	18 <sup>th</sup> Edition of Standard Methods 6410

**APPARATUS AND MATERIALS**

**Y      N**

1. Is glassware heated at 400°C (15-30 min.) or rinsed with pesticide grade acetone & hexane? [3.1] **NOTE:** Volumetric glassware should not be heated.
2. Are KD concentrator tubes calibrated (calibrations checked) at the volume/s employed? [5.2.3]
3. Have boiling chips been heated to 400°C for 30 min. for Soxhlet extract with methylene chloride? [5.3]
4. Is splitless injection employed for capillary columns? [5.6.1]
5. Is the data system capable of extracted ion current profile (EICP) plots--a plot of m/z vs. time? [5.6.6]


**SAMPLING AND STORAGE**

6. Are samples collected in glass containers (amber preferred) with lids lined with Teflon or foil? Automatic sampling equipment must be as free as possible of Tygon tubing. [9.1]
7. Are samples kept at 4°C from time of collection until extraction? [9.2]
8. Were samples dechlorinated using sodium thiosulfate at time of collection? [9.2]


**EXTRACTION - Separatory Funnel**

9. Is extraction completed within 7 days of collection and completely analyzed within 40 days of extraction? [9.3]
10. Is a reagent water blank processed with each set of samples or when reagents are changed? [8.1.3]
11. Is entire sample extracted? [10.2]
12. Is a surrogate standard spiking solution containing at least three surrogates at a concentration of 100µg/mL added to sep funnel prior to extraction with the pH adjusted to >11 using 10 N sodium hydroxide solution? [10.2]
13. What surrogates are being used?
14. Are samples extracted at pH >11 (NaOH) via three separate extractions--60 ml CH<sub>2</sub>Cl<sub>2</sub> each or extracted for 24 hours using continuous extraction? [10.2]
15. Are samples extracted at pH <2 (H<sub>2</sub>SO<sub>4</sub>) via three separate extractions--60 ml CH<sub>2</sub>Cl<sub>2</sub> each or extracted for 24 hours using continuous extraction? [10.5]
16. Is each fraction dried by passing it through NaSO<sub>4</sub>? [10.7]


Y N

17. Are the B/N and A extracts concentrated to about 0.5 ml using Kuderna Danish Apparatus (macro- then micro- Snyder columns)? [10.9]
18. Are the Snyder column balls wetted? [10.8]
19. Is the hot water bath temperature (60-65°C) and vertical position of the KDs adjusted to allow completion of concentration within 15-20 min. for 3-ball Snyder and 5-10 min. for 2-Ball Snyder column? [10.8 and 10.9]
20. Is final volume of extracts adjusted to 1.0 mL? [10.9]

#### **EXTRACTION - Continuous**

21. Is extraction completed within 7 days of collection and completely analyzed within 40 days of extraction? [9.3]
22. Is a reagent water blank processed with each set of samples or when reagents are changed? [8.1.3]
23. Is entire sample extracted? [10.2]
24. Is pH adjusted to >11 with 10 N sodium hydroxide solution and then transferred to continuous extractor? [11.2]
25. Is 1 mL of surrogate standard spiking solution added to extractor? [11.2]
26. Is sample bottle rinsed twice with 50-100 mL of methylene chloride? [11.2-3]
27. Is additional (200-500 mL) methylene chloride added to distilling flask along with sufficient water for proper operation and extracted for 24 hrs.? [11.4]
28. Is extract properly dried, concentrated and sealed? [11.4]
29. Is a clean distilling flask charged with 500 mL of methylene chloride, pH of aqueous phase adjusted to <2 (H<sub>2</sub>SO<sub>4</sub>), and extracted for 24 hrs.? [11.5]
30. Is extract properly dried, concentrated and sealed? [11.5]

#### **CALIBRATION**

31. Are stock standard solutions stored sealed in Teflon-sealed screw-cap bottles at 4°C (-10°C suggested)? [6.7.2]
32. Are stock standards replaced after 6 months? [6.7.3]
33. Does calibration involve a minimum of three conc. with one standard at the minimum reporting limit? [7.2.1]
34. Are internal standards used, does it contain at least 3 analytes from each group? [7.2]
35. Is the calibration curve verified each working day by analyzing at least one standard (near the expected sample conc.) with recovery of  $\pm 20\%$ ? [7.3]
36. Before analysis of samples, is a reagent water blank analyzed to demonstrate that system is under control? [8.1.3]
37. Are check standards analyzed at a rate of 5% of samples? [8.4]

#### **QA/QC**

38. Have the following operations been performed by each analyst to demonstrate accuracy and precision? [8.2]
  - a.) Were four quality control (QC) check samples containing each parameter of interest at a concentration of 100 µg/L in reagent water prepared and analyzed? These samples must have been obtained from a source separate from the calibration standards used. [8.2.2]

Y N

b.) Did the standard deviation (s) and average recovery (X) in µg/L meet acceptance criteria listed in Table 6? [8.2.5]

c.) If there was a failure, were the four aliquots prepared again and analyzed for the failed parameter? [8.2.6.2]

d.) If any parameter failed more than once, were all parameters reanalyzed? [8.2.6.2]

39. Are 5% of samples from each sampling site spiked? If only 1-20 samples are analyzed per month, only one spike is required. [8.3]
40. What is the concentration of the sample spike? (1 to 5 times background is recommended.) [8.3.1]
41. Do spike recoveries meet the acceptance criteria? [8.3.3]
42. For each failed spike, was a QC check standard analyzed and acceptable recovery (listed in Table 6) achieved? [8.4]
43. Is each sample, standard, and blank spiked with surrogate standard spiking solution (minimum of three surrogate compounds)? [8.6]
44. Are the DFTPP criteria in Table 9 met before any sample analysis is performed? [12.3]
45. Are the GC columns (packed) performance verified each day?  
B/N via Benzidine (tailing factor for 100 ng <3.0) [12.4]  
Acid via Pentachlorophenol (tailing factor for 50 ng <5.0) [12.5]
46. Is internal standard added to sample extract immediately before it is injected into the instrument? [13.3]
47. Is solvent flush technique used for manual injections, which are limited to 2-5 µL? (Smaller volumes (1.0 µL) may be injected using an autosampler) [13.4]
48. Does the retention time of the sample target compound agree within  $\pm$  30 seconds of that measured for the standards? [14.1.2]
49. Do the relative intensities of the three characteristic ions agree within  $\pm$  20% of the relative intensities listed in Table 4 & 5? [14.1.3]
50. Are the compound concentrations in the sample calculated using the response factor (RF)? [15.1]  
NOTE: If the RF is <35% RSD, the average RF may be used. [7.2.1]
51. Are samples diluted and re-analyzed if the measured conc. is above the highest calibration std. conc.? [13.5]

PROBLEMS: